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Material properties

**Ref :
282 027**

Français – p 1

English – p 5

Version : 9001

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Immersible resistances

1 Principe – description

1.1 But

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1 General description

1.1 Purpose

The device is mainly intended to be immersed in a calorimeter. It allows for:

- demonstration of the formula giving the amount of heat dissipated by a current flow through a resistance.
- determination of the value of J.
- measurement of the specific heat of various liquids.

1.2 Characteristics

The device can be used with any calorimeter but has been especially designed to fit to the simple calorimeter Ref. 251 014 supplied with its cover Ref. 251 015. It contains 2 resistances of 2 and 4 Ω coiled with Nichrome wire \varnothing 3/10 (80% Ni; 20% Cr; resistance 15.3 Ω /m), adjusted at $\pm 10\%$. The 2 resistances are connected in series in such a way that by using the 3 output terminals as is indicated on the support, it is possible to obtain the 3 values: 2 Ω , 4 Ω , 6 Ω .

Temperature of the new wire \varnothing 0.3 mm stretched in quiet air at 20°C					
1 ampere	:	200 °C	3 amperes	:	850 °C
2 amperes	:	550 °C	4 amperes	:	1100 °C

1.3 Operation

Water can be used, preferably deionised, but the use of petrol is recommended for the majority of the experiments: low specific heat, no electrolysis.

The device is dipped in the calorimeter; the top flange bearing the terminals completely blocks the calorimeter orifice off. The width of the rod base is such that the stirrers that equip this calorimeter cannot be in contact with these rods.

The device is handled by holding it by the top flange, or even more simply by using inlet plugs of the current as the load-handling device.

2 Experiments

❖ Material required for the experiments

- 282 027 Immersible resistances.
- 351 036 Digital chronometer.
- 281 009 Rheostat 33 Ω 3.1 A.
- 291 133 Multimeter.
- 251 014 Simple calorimeter.
- 251 015 Cover with stirrer.

- 281 083 Power supply 6 V – 12 V
- 253 067 Thermometer -10 + 60 °C.

2.1 Variation of the amount of heat with the amperage

The water equivalent m of the system (calorimeter, thermometer, stirrer and resistance) needs only be determined once.

Comments:

The aluminium alloy rods weigh 12 g, corresponding to an m value of around 5.

The nuts are made of aluminium to avoid corrosion by formation of hydro electrical coupling.

Place a weight of petrol M (of around 300 g) whose water equivalence is MC (C being the specific heat of petrol) in the cell of the calorimeter. Immerse the resistance and the accessories of the calorimeter, and record the temperature of the unit displayed on the thermometer. Constitute a series circuit containing: a current source, an ammeter, a rheostat and the resistance.

Close the circuit and begin timing. Record the amperage displayed on the ammeter. Let the current flow during a time T . Run the experiment again several times with amperage modifications.

Amounts of heat: $Q = (MC + m) (\theta_2 - \theta_1)$

M is the weight of the petrol.

C is the specific heat of the petrol.

M is the water equivalence of the calorimeter and its accessories.

θ_1 is the initial temperature.

θ_2 is the final temperature.

If I_1 , I_2 and I_3 are the different measured amperages, it must be checked that the Q_1 , Q_2 ...etc. ratios of the emanated amount of heat are between them as should the squares of the measured amperages be.

2.2 Variation of the amount of heat with time

Record the temperature θ_1 several times at frequent intervals. Calculate the amount of heat emanated: $Q = (MC + m) (\theta_i - \theta_1)$

It is observed that the amounts of heat emanated are proportional to the times separating the beginning of the experiment from each of the readings.

2.3 Variation of the amount of heat with resistance

Run the experiment again with the same amperage and time conditions, successively using the following resistance values: R_1 , R_2 and R_3 .

2.4 Determination of the W ration (J calculation)

Let the current flow during a time of T seconds. The amount of energy emanated is $W = RI^2T$ in Joules. By calculating the W/Q ratio, you must find a number approximately 4 ($J = 4.18$)

2.5 Measurement of the specific heat of liquids

One interesting application of Joule's law is to measure the specific heat of various liquids (preferably organic ones, to avoid the attack of the metallic parts).

A temperature increase of θ_1 is obtained. Replace the water by another liquid and run the experiment again with the same resistance, the same amperage and the same time. The temperature increase is θ_2 . The liquid specific heat M is calculated with the formula:

$$C/1 = \theta_1 / \theta_2$$

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries